

Article-Number: 8053 100 002

Typ: XELSIUS Reactor

Manufacturer:

nevoLAB GmbH  
Am Gehrenbach 8  
D 88167 Maierhöfen  
www.nevolab.de

Telefon:  
+49 8383 9929 566 0

# XELSIUS Cooling Guide

## Understanding the XELSIUS System

| Version | Date       | Subject | Abr. |
|---------|------------|---------|------|
| 1.1     | 14.01.2020 | Review  | MJ   |
| 1.0     | 28.11.2019 | Review  | ML   |
| 0.1     | 25.11.2019 | Initial | ML   |

|  |   |
|--|---|
| <i>Content</i>                               | 1 |
| A Preface                                    | 2 |
| B Safety Information                         | 3 |
| <i>Safety Instructions</i>                   | 3 |
| C Learn about XELSIUS and the cooling system | 4 |
| <i>XELSIUS</i>                               | 4 |
| <i>Factors defined by Reaction</i>           | 5 |
| <i>Cooling Liquid</i>                        | 5 |
| <i>Parameter Summary</i>                     | 5 |
| <i>Working Temperature Range</i>             | 5 |
| <i>Cooling Capacity</i>                      | 5 |
| <i>Flow Rate</i>                             | 5 |
| <i>Pressure</i>                              | 5 |
| <i>Other</i>                                 | 5 |

## Content

## A Preface

This manual helps to commission and to operate the device safely and smoothly. The instructions must therefore be read carefully before use. Read this manual before using the equipment.

Please keep this manual accessible to all users at all times.

In these operating instructions, the following safety symbols are used. They are intended to alert the reader to the associated safety instructions in the text.



This symbol indicates risks to people, equipment, materials or the environment

The documentation has been prepared with care. The product described is constantly evolving. Therefore, the documentation is not in every case checked for consistency with performance data, standards or other characteristics are not a form of guarantee in legal aspects.

If it contains technical or editorial errors, we reserve the right to make alterations at any time and without notice.

© This manual is copyrighted. Any reproduction or third party use of this publication, in whole or in part, is prohibited without written permission of the manufacturer.

## B Safety Information



- The XELSIUS Reactor must be used with an appropriately rated cooling system at all times, even when the equipment is only used for heating applications.
- Failure to use an appropriately rated cooling system will negatively impact the performance of XELSIUS. It could reduce the lifespan of the product, increase service frequency or cause damage of the equipment.
- Damage caused to the XELSIUS Reactor by using inappropriately rated cooling systems is not covered under warranty.
- Damage caused to the XELSIUS Reactor by using inappropriately cooling liquides (unfiltered, contaminated or wrong liquids) is not covered under warranty.

## Safety Instructions



## C Learn about XELSIUS and the cooling system

XELSIUS is capable of independently cooling, heating and stirring each of its 10 cells down to a minimum temperature of at least  $-20\text{ }^{\circ}\text{C}$ , but this is only achievable when the product is coupled to a cooling system that can remove heat at a rate greater or equal to that collected by the XELSIUS system.

It is important to understand the importance of this relationship and to select an appropriate cooling system so that you can achieve the performance needed for your application and get the results you expect.

Each of the XELSIUS Reactor cells contains its own Peltier cooling/heating device. While operating in cooling mode, these devices remove heat from the reaction vessel and pass it onto the cooling system via a heat exchanger.

The temperature of the coolant increases which is then pumped out of the reactor, where it's cooled by the chiller/cooling system and rejected into the air.

It is important to note that even when a cell is heating and not using the cooling feature, a cooling system is still required to ensure damage is not sustained to the sensitive cells during operation.

If the cooling system is unable to transport the heat produced inside the system in time, the performance (minimum temperature/ cooling rates) will be effected negatively.

XELSIUS

Factors which contribute to the load on the cooling system:

- The number of cells being cooled  
Cooling 10 cells in parallel places a higher load on the cooling system than only cooling a lower number in parallel
- The cooling rate  
Cooling at a rate of 5 °C/minute places a higher load on the cooling system than cooling at 1 °C/minute.
- The minimum temperature required.  
Cooling down to a temperature of -25 °C places a higher load on the cooling system than cooling down to a temperature of 0 °C
- Heat of the reaction  
Strong exotherm reactions places a higher load than endotherm reactions

Once you have selected your cooling system, you will need to follow the manufacturer's recommendation and use a fluid suitable for your application. Only use anti-freeze fluids designed specifically for lab equipment.

When working with refrigerated and other cooling equipment, it is important to use a fluid that will not freeze because the temperature inside the cooling equipment may be significantly colder than the set temperature.

Water-based synthetic fluids are popular choices for cooling. They can be used across a broad temperature range. Glycols, when mixed with water, become antifreeze. There are two glycols which are commonly used: Ethylene glycol (EG) and propylene glycol (PG). Of the two, EG has slightly better thermal properties and can cool in a broader temperature range. However, it is highly toxic and not usable in all applications. In those cases, PG is a good alternative. We recommend mixing with distilled water.

Factors defined  
by Reaction

Cooling Liquid



The main parameters for a cooling system are:

- Working Temperature Range  
typically: 1 °C – 30 °C
- Cooling Capacity Chiller  
typically: 1 kW in the range of 0 °C – 15 °C
- Flow Rate  
minimum flow rate 3,0 l/min  
maximum flowrate 20,0 l/min
- Pressure  
specified between 1,0 bar and 2,0 bar
- Filtering Requirements  
Inline filter: 50 micron
- Cooling Liquid  
water or water/glycol mixtures

The main parameter to consider when defining a cooling system is the working temperature range.

The working temperature you use will be subjective to your experiment and is typically determined empirically.

Here are some parameter settings from our practical experiences with the XELSIUS system.

Chiller: Huber, Unichiller 010 Olé



Flow rate: 10,0 l/min

Ambient temperature: 23 °C

Running 10 cells parallel with max cooling / heating temperature.

Parameter  
Summary

Working  
Temperature  
Range



**Results**

Maximum Cell Temperatures

| Chiller Temperature [°C] | Cell [°C] | +/- | Criterion |
|--------------------------|-----------|-----|-----------|
| 35,0                     | 155,0     | 1,5 | > 150     |
| 20,0                     | 141,0     | 1,0 | > 140     |
| 8,0                      | 134,0     | 1,5 | > 120     |
| 1,0                      | 126,0     | 1,5 |           |

Minimal Cell Temperatures

| Chiller Temperature [°C] | Cell [°C] | +/- | Criterion |
|--------------------------|-----------|-----|-----------|
| 40,0                     | 4,7       | 0,8 | < 0       |
| 20,0                     | -9,6      | 1,0 |           |
| 13,3                     | -17,6     | 1,5 | < -15,0   |
| 12,0                     | -19,0     | 1,6 |           |
| 10,0                     | -20,5     | 1,5 |           |
| 8,0                      | -21,5     | 2,5 | < -20,0   |

How to use the diagram?

With a working temperature of 20 °C, (here provided by a Huber Unichiller 010 Olè with a minimum flowrate of 3,0 l/min ) and running the XELSIUS with 10 cell parallel the reachable jacket temperatures will be:

lowest temperature: -10 °C  
highest temperature 140 °C

With a working temperature of 8 °C cooling of the jacket will be at approx. -21 °C, while maximum temperature will be approx 134 °C.

The cooling capacity is a characteristic of a chiller which determines the performance and the speed of the cooler. The cooler needs a capacity greater or equal than that produced within the reactor.

It is recommended that you select a product with a cooling capacity of at least 1.0kW (1000W) between the working temperature of range of 15 °C to - 5 °C

Cooling Capacity: for example, Huber Unichiller 010 Olé

|                  |        |        |        |         |
|------------------|--------|--------|--------|---------|
| Working Temp     | 15 °C  | 0 °C   | -10 °C | -20 °C  |
| Cooling Capacity | 1,0 kW | 0,8 kW | 0,5 kW | 0,15 kW |

It is recommended to use the warmest working temperature possible which achieves your performance requirements to avoid unnecessary condensation and to avoid wasting energy.

If the cooling capacity of your cooling system is not strong enough your cooler will not be able to hold the set temperature.

Please note that the cooling capacity decreases significantly when the working temperature is reduced.

Flow Rate is also a very important parameter for your cooling system. Physically determined heat exchange is like other transport processes a kinetic phenomenon. The cooler should be capable of providing coolant at a rate of

typically minimum flow rate 3,0 l/min and maximum flowrate 20,0 l/min

Higher flow rates are useful to reach lower jacket temperatures, while lower flow rates are favorable for higher temperatures.

In every case you must ensure that the supply of 3 Litres per minute is maintained for the duration of your work, even if you share your cooling system with other devices.

Do not compromise the flow rate of the coolant into the XELSIUS by using incorrectly specified ancillary items such as couplers, reducers or tee-pieces, for example.

## Cooling Capacity

## Flow Rate



## Learn about XELSIUS and the cooling system

In order to obtain the flow rate discussed above, the coolers pump must be capable of providing a pressure of at least 1 bar effective at the XELSIUS system.

You must ensure that the pressure of at least 2 bar is maintained for the duration of your work.

If the flowrate and pressure are not effective to the XELSIUS system it may cause overtemperatures within the reactor system.

At a temperature higher than 60 °C inside the reactor, due to safety reasons the power supply for the cells will be terminated.

If you use colling system with lower specifications it may be that you could not reach the full potential of the XELSIUS and may not meet the performance characteristics required by your work.

When coupling your XELSIUS to a a cooling system it is required that you also fit a 50 micron inline filter (if one is not present in the cooler) to prevent the ingress of contaminates into the cell heat exchangers. The filter should be installed on the outlet of the chiller.

Partial restriction or blockage of heat exchangers caused by contaminates will impact the system performance and potentially void your warranty.

### Pressure



### Other

