


Data Devices



nevoLAB Link

Art.: 8027 420 001 1 pcs.

- Connect the robot to your computer wirelessly
- Range of ca. 30-100m
- 1m USB-cable



nevoLAB RFID

Art.: 8027 440 001 1 pcs.

127 x 85 x 17mm, white

- Read and write RFID tags on the vials
- Fits into a rack position on the Samplify Robot
- 1m USB-cable



RFID-Pen

USB

Art.: 8027 440 002 1 pcs.

Wireless

Art.: 8027 450 001 1 pcs.

- Read and write RFID tags on vials or racks
- Works with ultra-small RFID-tags



Barcode Reader

Art.: 8027 460 001 1 pcs.

- Read out barcodes on racks and vials
- 1D & 2D-barcodes
- 1m USB-cable

Key Features

- Works with all SAMPLIFY devices
- Can be easily controlled in Samplify Script



nevoLAB GmbH
Am Gehrenbach 8
D 88167 Maierhöfen
phone: +49 (0) 8383 929 566 0
fax: +49 (0) 8383 929 566 30
email: info@nevolab.de
www@nevolab.de

Version: 2016-08-24