

Racks for HPLC/GC vials

SAMPLIFY can be equipped with various racks. Following labware is tested and available in our SAMPLIFY Script database.



Rack 6 x 9

Vial plate for 54 x 2 mL vials
Art.: 8027 330 001 1 pcs.



Rack 3 x 5

Vial plate for 15 x 6 mL vials
Art.: 8027 330 002 1 pcs.



Rack 3 x 9

Vial plate for 27 x 2 mL Eppendorf tubes
Art.: 8027 330 003 1 pcs.



Rack 6 x 8

Vial plate for 48 x 2 mL vials, Waters
Art.: 8027 330 011 1 pcs.



Rack 3 x 4

Vial plate for 12 x 6 mL RFID-vials
Art.: 8027 330 012 1 pcs.




Rack 4 x 5 Shaker

Vial plate for 20 tubes with Ø 19 mm, aluminium, usable with shaker
Art.: 8027 330 018 1 pcs.


Key Features

Vial plates
for volumes up to 2 or 6 mL
for HPLC and RFID-vials



nevoLAB GmbH
Am Gehrenbach 8
D 88167 Maierhöfen
phone: +49 (0) 8383 929 566 0
fax: +49 (0) 8383 929 566 30
email: info@nevolab.de
www@nevolab.de

Version: 2016-08-24



ANSI/SBS Well Plates

SAMPLIFY can be equipped with various racks. Following labware is tested and available in our SAMPLIFY Script database.



MTP 96

96-Well Assay Microplate, PP, ANSI/SBS 1-2004
Art.: 8027 330 004 1 pcs.



24 DeepWell

24-DeepWell Assay Microplate, PP, for Ø 16 mm tubes
ANSI/SBS 1-2004
Art.: 8027 330 014 1 pcs.



96 DeepWell

96-DeepWell Assay Microplate, PP, ANSI/SBS 1-2004
Art.: 8027 330 015 1 pcs.



6 Well Cell Culture Plate

6-Well Cell Culture Plate, PS, ANSI/SBS 1-2004
Art.: 8027 330 016 1 pcs.




24 Well Cell Culture Plate

24-Well Cell Culture Plate, PS, ANSI/SBS 1-2004
Art.: 8027 330 0171 pcs.


Key Features

Microplates and
Cell Culture Plates



nevoLAB GmbH
Am Gehrenbach 8
D 88167 Maierhöfen
phone: +49 (0) 8383 929 566 0
fax: +49 (0) 8383 929 566 30
email: info@nevolab.de
www@nevolab.de

Version: 2016-08-24



Tube Racks

SAMPLIFY can be equipped with various racks. Following labware is tested and available in our SAMPLIFY Script database.




Test Tube Racks

Test tube racks, PP, autoclavable, 265 x 126 x 75 mm, usable with Adapter 8027 310 006

Art.: 8027 330 005	3x7 Tubes with Ø 30 mm	1 pcs.
Art.: 8027 330 006	4x8 Tubes with Ø 25 mm	1 pcs.
Art.: 8027 330 007	4x10 Tubes with Ø 20 mm	1 pcs.
Art.: 8027 330 008	5x11 Tubes with Ø 16 mm	1 pcs.
Art.: 8027 330 009	5x11 Tubes with Ø 18 mm	1 pcs.
Art.: 8027 330 010	6x14 Tubes with Ø 13 mm	1 pcs.

Key Features

Test Tube Racks
for volumes up to 50 ml



Drip Tray

Drip Tray, HD-PE, 127 x 85 x 18 mm
Art.: 8027 330 013 1 pcs.



Rack 2 x 4

Vial plate for 8 tubes, Ø 28 mm, PP, black, ANSI/SBS 1-2004
Art.: 8027 330 019 1 pcs.



nevoLAB GmbH
Am Gehrenbach 8
D 88167 Maierhöfen
phone: +49 (0) 8383 929 566 0
fax: +49 (0) 8383 929 566 30
email: info@nevolab.de
www@nevolab.de

Version: 2016-08-24