



## Needle Wash Units



### Key Features

- Compatible with all SAMPLIFY Robots
- Chemically resistant PTFE

ID	Qty	Description
8027 240 001	pcs	Washstation: PTFE, ID 5,2 mm, optimized for 1/8"-tubes Hosebarb: PP, 10mm
8027 240 002	pcs	Washstation: PTFE, ID 3 mm, optimized for needles with 1,4 mm diameter Hosebarb: PP, 2,5mm
8027 240 003	pcs	Washstation: PTFE, ID 2 mm, optimized for needles with 1,4 mm diameter Hosebarb: PP, 2,5mm



nevoLAB GmbH  
Am Gehrenbach 8  
D 88167 Maierhöfen  
phone: +49 (0) 8383 929 566 0  
fax: +49 (0) 8383 929 566 30  
email: [info@nevolab.de](mailto:info@nevolab.de)  
[www@nevolab.de](http://www@nevolab.de)

Version: 2016-08-24