

# Starter Kits & Essentials



**xelsius**


TEN REACTIONS ⇌ ONE DEVICE



## Reflux Condenser Modul

Item No 8053 000 201

- (W/H/T) 352 mm / 156 mm / 48 mm
- Material: PSI
- Material cooling body: Stainless steel
- Temperature range: 5 °C to 35 °C
- Flow rate: 1,2 l/min to 5 l/min
- Pressure: max. 2 bar




## HV Vials Starter Kit | 1 – 30 ml

Item No 8053 000 202

- GL 25 Standard
- Volume: 1–30 ml

10 × HV Vials GL25  
 10 × Magnetic stirring bar  
 10 × Septa N24  
 10 × Cap GL25




## LV Vials Starter Kit | 0,5 – 5 ml

Item No 8053 000 203

- GL 25 Standard
- Volume: 0,5–5 ml

10 × LV Vials GL25  
 10 × small Magnetic stirring bar  
 10 × Septa N24  
 10 × Cap GL25  
 10 × Adapterset reactor for LV Vials



## On Top Adapter xelsius

Item No 8053 800 023

- Material: PE
- Compact build up for xelsius
- Compatible for xelsius

All information without any warranty