

# Temperature Measurement



**xelsius**

TEN REACTIONS ⇌ ONE DEVICE



## Liquid Temperature Probe | Set of 5

Item No 8053 000 204

- Temperature probe PT 1000
- Temp. range: -30 °C to 170 °C
- Accuracy: ±0,8 °C
- 5 × Liquid temperature probes, ready to connect into xelsius.
- 5 × PT 1000 with Lemo Connector, FEP coated (-50 °C to 200°C)
- 5 × PTFE adapter and UNF 28 1/4" fitting
- 10 × Pre-punched GL25 septa



All information without any warranty