



xelsius

TEN REACTIONS ⇌ ONE DEVICE


# 5-Port-Reaction Cap




## 5-Port-Reaction Cap and Vialset

Item No 8053 000 205

- GL 25 standard
- Volume: 1–30 ml
- 1 × HV Vials GL 25 1–30 ml
- 1 × Magnetic stirring bar
- 1 × 3,2 mm fitting
- 1 × 1,6 mm fitting
- 4 × Blind plug
- 2 × Tube clip 2,4 mm




All information without any warranty



### 5-Port-Reaction Cap Config Set

Item No 8053 000 207

- 5 × 5-Port-Reaction Cap
- 10 × 3,2 mm fitting
- 5 × 1,6 mm fitting
- 20 × Blind plug
- 2 × Tube clip 2,4 mm
- 10 × Luer male
- 10 × Luer female tube clip
- 10 × Septa N9, PTFE / silicone
- 100 × Septa N9 preslit, PTFE / silicone
- 2 × Bidirectional-valve



### Needle Research Kit

Item No 8053 000 209

- for multiple applications in sampling and dosing together with 5-Port-Reaction Cap
- 2 × Disposal needles  
0,8 × 120 mm 21 G × 4 3/4"
- 1 × Universal needle SS18 6" 150 mm
- 1 × Stainlesssteel needle PNC20 6", 150 mm
- 1 × FEP tube 3,2 mm w/ Luer Tip and PEEK Fitting 3,2mm
- 1 × FEP tube 1,6 mm w/ Luer-Tip and PEEK Fitting 1,6 mm
- 1 × Luer syringe 5,0 ml
- 1 × Luer syringe 1,0 ml
- 100 × N9 Septum Silicon preslit



### Set of 10 GL25 Caps

Item No 8053 800 027

- GL25 cap with hole
- Material: PS
- with septa PTFE/silicone 3,2 mm

All information without any warranty