


Inerting Workstation & Gas Injecting



xelsius

TEN REACTIONS ⇌ ONE DEVICE



Inert Workstation for 5-Port-Reaction Cap

Item No 8053 000 206

- Compatible for nitrogen & inert gas
- Pressure: -1 to 1 bar
- Material: Stainless steel

10 × HV Vials; GL 25, 1–30 ml

10 × 5-Port-Reaction Cap

10 × Magnetic stirring bar

20 × Blind plug

20 × Tube clip 2,4 mm

1 × 1 m tube silicone,
Ø 6/2 mm



Gas-Balloon Supply Set

Item No 8053 000 208

- for easy and safe gas injecting

1 × Ballonadapter

1 × 2/2 PTFE Valve

1 × Tubing kit with selfsealing couplings

1 × 5-Port-Reaction Cap w/ fittings

2 × PEEK fitting 1,6 mm w/ FEP Tube

10 × Ballons



Tube Coupling, Ø 3,2 mm | Set of 5

Item No 8053 800 028

- Self sealing
- Material: PP
- Seal: NBR
- Tube port: ID = 3,2 mm
- Non spile male and female

All information without any warranty