

Reaction Vessels & Adapters



HV Vial, GL 25, 1–30 ml

Set of 10

8053 800 004

134 × Ø 24mm, GL 25 Standard,
Volume: 1–30 ml



LV Vial, GL 25, 0,5–5 ml

Set of 10

8053 800 005

134mm × Ø 24 mm | 13 mm, GL 25 Standard
Volume: 0,5–5 ml



HV Vials Water Separator

Set of 10

8053 800 041

150 × Ø 24 mm, GL 25 Standard,
Volume: 1–30 ml, Separator: 2 ml



HV Vial PTFE, GL 32

Set of 5

8053 800 021

GL 32 Standard, PTFE, Volume: 1–30 ml



5-Port-Reaction Cap and Vialset

8053 000 205

1 × HV Vials GL 25 1–30 ml
1 × Magnetic stirring bar
1 × 3,2 mm fitting
1 × 1,6 mm fitting
4 × Blind plug
2 × Tube clip 2,4 mm



Multiblock NMR Tube

Set of 4

8053 800 075

5 Position per block, anodized aluminum
for NMR tubes; OD 5 mm



Multiblock Vial OD 8 mm

Set of 4

8053 800 079

4 Position per block, anodized aluminum
for vials; OD 8 mm



Adapter-Vials 24/19 mm

Set of 10

8053 800 046

Adapter to use individual vessels
with OD: 19 mm, other OD on request