

# Certificate

## 8053 000 211 - xelsius Pressure Reactor 15 bar

Articel number	8053 000 211 / xelsius Pressure Reactor 15 bar
Serial number	53-211-001-23
Manufacturer	nevoLAB GmbH Am Gehrenbach 8 D-88167 Maierhöfen
Applied directive	Directive 2014/68/EU



### Checklist

Connection between pressure cap - Pressure Monitoring & Safety Modul	<input checked="" type="checkbox"/>
Assambly safety Ventil & pressure gauge	<input checked="" type="checkbox"/>
Assambly Inlay PTFE 10 ml	<input checked="" type="checkbox"/>
Assambly Inlay Glass 7 ml	<input checked="" type="checkbox"/>

### Pressure stability test

Test medium	Water H2O
Test pressure	40 bar
The pressure reactor has passed the stability test <input checked="" type="checkbox"/>	

### Pressure leak test

Test medium	nitrogen N2
Test pressure	20 bar
Pressure drop after 1h	max. 0,1 bar
The pressure Reactor has passed the leak test <input checked="" type="checkbox"/>	

### Comment

Maierhöfen, Date:	Inspector	Signature
27.06.2023	HB	Brand